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¹Frans van Wagenburg et al. UT Southwestern. TAVR Handbook: Step-by-step guide: Transcatheter aortic valve replacement. American College of Cardiology, <https://www.acc.org>. 2018; Chapter 15

²Chamandi, C et al. Transcatheter Aortic Valve Replacement Compared With Other Alternative Access Routes For Transcatheter Aortic Valve Replacement. 16 Nov 2018 <https://doi.org/10.1161/CIRCINTERVENTIONS.118.006388>
Circulation: Cardiovascular Interventions. 2018;11:e006388

³Kirker, Eric B. et al. The Carotid Artery as Preferred Alternative Access Route for Transcatheter Aortic Valve Replacement. The Society of Thoracic Surgeons, published by Elsevier Inc. 2107.
[https://www.annalsthoracicsurgery.org/article/S0003-4975\(16\)31888-4/pdf](https://www.annalsthoracicsurgery.org/article/S0003-4975(16)31888-4/pdf)

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- Outcomes Following Transcarotid Versus Transapical and Transaortic Access for TAVR with Balloon-Expandable SAPIEN 3 Valve - An analysis of the STS/ACC Transcatheter Valve Therapy (TVT) Registry Keith B. Allen¹, Adnan K. Chhatriwalla¹, John Saxon¹, James Hermiller², David Heimonsohn², Sina Moainie², Raymond Mckay³, Mohiuddin Cheema³, Brandon Jones⁴, Robert Hodson⁴, Ethan Korngold⁴, Eric Kirker⁴ Presented at TVT 2019
- Trans-carotid Versus Subclavian/Axillary Access for Transcatheter Aortic Valve Replacement: Real World Comparative Outcomes Utilizing the STS/ACC TVT Registry James Hermiller¹, David Heimansohn¹, SinaMoainie¹, Eric Kirker², Ethan Korngold², Robert Hodson², Brandon Jones², Adnan K. Chhatriwalla³, John Saxon³, Keith B. Allen³, Raymond Mckay⁴, Mohiuddin Cheema⁴ Presented at TCT 2019
- Transcarotid Versus Subclavian/Axillary Access for Transcatheter Aortic Valve Replacement: A Propensity-Matched Analysis of Real World Data Eric Kirker¹, Ethan Korngold¹, Robert Hodson¹, Brandon Jones¹, Raymond Mckay², Mohiuddin Cheema², David Heimansohn³, SinaMoainie³, James Hermiller³, Adnan K. Chatriwalla⁴, John Saxon⁴, Keith B. Allen⁴ Presented at STS 2020
- Frans van Wagenburg et al. UT Southwestern. TAVR Handbook: Step-by-step guide: Transcarotid, transcaval and direct aortic. American College of Cardiology, <https://www.acc.org>. 2018; Chapter 15
- Chamandi, C et al. Transcarotid Compared With Other Alternative Access Routes For Transcatheter Aortic Valve Replacement. 16 Nov 2018 <https://doi.org/10.1161/CIRCINTERVENTIONS.118.006388> Circulation: Cardiovascular Interventions. 2018;11:e006388
- Kirker, Eric B. et al. The Carotid Artery as Preferred Alternative Access Route for Transcatheter Aortic Valve Replacement. The Society of Thoracic Surgeons, published by Elsevier Inc. 2017. [https://www.annalsthoracicsurgery.org/article/S0003-4975\(16\)31888-4/pdf](https://www.annalsthoracicsurgery.org/article/S0003-4975(16)31888-4/pdf)

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©2024 Edwards Lifesciences Japan LLC. All rights reserved. 無断転載を禁じます EW2024065

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