


# The ReTAVI Trial: Clinical Outcomes of Patients Undergoing a Redo-TAVI Procedure

G. Tarantini, R. Parma - Evaluation of Clinical Outcomes of Patients Undergoing a Redo-TAVI Procedure; a Multicenter Prospective Observational Registry. Presented at TCT 25<sup>th</sup> October 2025, San Francisco, USA.



| Study design                                                                            | Device(s) used                                                                            |                                                                                                                                                       | Locations                                               | Patients characteristics                                         |
|-----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|------------------------------------------------------------------|
| <br>Prospective Observational Registry<br><b>N=150</b><br>Follow up:<br>30d, 1y, 3y, 5y | Index valve<br><br>The SAPIEN platform<br><b>30.0%</b>                                    | Second valve<br><br>SAPIEN 3<br>SAPIEN 3 Ultra<br>SAPIEN 3 Ultra RESILIA<br>Sizes used:<br>20 mm: 8.4%<br>23 mm: 39.2%<br>26 mm: 44.8%<br>29 mm: 7.7% | <br><b>59 sites</b><br>Across Europe, Canada and Israel | <br>Median Age: <b>84 years</b><br>Median STS score: <b>7.0%</b> |
|                                                                                         | <br>Core valve/ Evolut<br><b>53.1%</b><br>ACURATE<br><b>14.0%</b><br>Other<br><b>2.8%</b> |                                                                                                                                                       |                                                         |                                                                  |

First prospective, multicenter international registry with patients treated with the SAPIEN 3 platform for THV-in-THV procedures



Proven excellent performance of the SAPIEN 3 platform in THV-in-THV procedures




Excellent outcomes for THV-in-THV procedure at discharge and 30 days

Key patient outcomes included (n=143):


|  |                               |                          |                        |
|--|-------------------------------|--------------------------|------------------------|
|  | <b>VARC-3 Device success</b>  | <b>95.0%</b>             |                        |
|  | <b>All-cause mortality</b>    | Discharge<br><b>3.5%</b> | 30 days<br><b>3.5%</b> |
|  | <b>All Stroke</b>             | Discharge<br><b>0.0%</b> | 30 days<br><b>0.7%</b> |
|  | <b>Pacemaker Implantation</b> | Discharge<br><b>5.6%</b> | 30 days<br><b>6.3%</b> |

Significant patient functional status improvement

NYHA Class scores at 30 days (n=143):



## 98.6% rate of no coronary obstruction with the SAPIEN 3 platform in studied index valves




|                            | Short index frame<br>(n=46) | Tall index frame<br>(n=97) |
|----------------------------|-----------------------------|----------------------------|
| Use of coronary protection | 22.2%                       | 28.3%                      |

Coronary protection  
was more used  
in tall frame valves

## Consistent hemodynamic improvement across anatomies, index valve types and sizes



Hemodynamic improvements in all SAPIEN 3 platform sizes



## The SAPIEN 3 platform demonstrated outstanding regurgitation results



Demonstrating absence of Intra-prosthetic regurgitation (n=143)

| Outcomes                                                | 30 days |
|---------------------------------------------------------|---------|
| Intraprosthetic regurgitation; <i>moderate/severe</i>   | 0.0 %   |
| Paravalvular regurgitation/leak; <i>moderate/severe</i> | 0.9 %   |

The image of the SAPIEN 3 Ultra RESILIA valve in the SAPIEN frame is for illustrative purposes only and does not reflect the study design, which includes other index valves. SAPIEN is an older model that has been discontinued and it is not available for commercial use

**Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult [eifu.edwards.com](http://eifu.edwards.com) where applicable).**

Edwards, Edwards Lifesciences, the stylized E logo, RESILIA, SAPIEN, SAPIEN XT, SAPIEN 3, and SAPIEN 3 Ultra are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2025 Edwards Lifesciences Corporation. All rights reserved. PP--EU-11334 v1.0

Edwards Lifesciences Sàrl • Route de l'Etraz 70, 1260 Nyon, Switzerland • [edwards.com](http://edwards.com)



Edwards