

CT screening measurements quick reference

CT imaging and analysis for transcatheter heart valve sizing

Overview

After locating the annular plane, measurements of the annulus and related structures are obtained. These measurements can be grouped into those that are dependent on the annular plane and those that are not.



Measurements dependent on the annular plane

- Annular dimensions
- Distance to the coronary ostia (Left and right coronary artery height)
- Sinus height (sinotubular junction heights)
- Ascending aorta length
- Left ventricular outflow tract diameter
- Prediction of fluoroscopic angle

Measurements independent of annular plane

- Sinus of valsalva width
- Sinotubular junction (STJ) diameter
- Ascending aorta diameter

Note: Although several automated tools are commercially available to expedite the measurement work flow, it is important to understand the approach for each measurement in order to be capable of performing the assessment manually, and to recognize errors which may occur in an automated analysis.

CT report quick reference

Reporting

The analysis report should contain the following data:

Acquisition

- Data acquisition mode
- Timing of images in the cardiac cycle (systolic vs diastolic)
- Contrast volume
- Image quality

Aorta

- Presence of kinking
- Presence of intraluminal obstruction
- Presence of intraluminal thrombus

Ascending aorta

- Width at 40 mm from annulus
- Position relative to sternum

Aortic arch

- Width
- Branch anatomy (for embolic protection device purposes and/or alternative access)

Descending aorta

- Width

Iliofemoral arteries

- Minimum width on both sides
- Tortuosity
- Calcification

Aortic root

- Sinotubular junction diameter
- Sinus of valsalva width
- Sinus of valsalva height
- Distance of coronary ostia from aortic annular plane

Aortic valve

- Cuspidity
- Qualitative extent of aortic valve calcification, separately for commissures and annulus
- Presence of a severely calcified cusp which may obstruct a coronary ostium

Aortic annulus

- Long and short diameters
- Area and area-derived diameter
- Perimeter and perimeter-derived diameter

Fluoroscopic coplanar angle

- C-arm angulation with right coronary cusp in the middle

Left ventricle

- Presence of thrombus
- Septal hypertrophy

Timing of images

It is important to report the phase in which the measurements were obtained due to the dynamism of the aortic annulus.

Image quality

It is important to report if there were acquisition and/or reconstruction artifacts in the dataset that could affect the accuracy of the provided measurements.

Aorta

If there is a significant amount of aortic disease, an alternative access route may be considered.

Iliofemoral arteries

If there is a significant amount of iliofemoral disease an alternative access route may be considered. Additional focus should be placed on the access site for the presence of calcification.

Aortic root, aortic valve and aortic annulus

All of the items listed in these categories are critical in the determination of the appropriate prosthesis size and optimal positioning of the prosthesis in the native anatomy, as well as the suitability of the patient for transcatheter aortic valve placement.

Fluoroscopic coplanar angle

The reporting of the predicted fluoroscopic coplanar angle can significantly decrease the amount of contrast required during the procedure to identify the deployment angle.

Left Ventricle

The presence of significant septal hypertrophy may alter the desired landing zone for the prosthesis.

This reference card is provided as an educational resource to medical personnel by Edwards Lifesciences (the "Author"). The information in this reference card has been compiled from then-currently available literature. Although every effort has been made to faithfully report the information and keep it up to date, the Author cannot be held responsible for the completeness or accuracy. This reference card is not intended to be and should not be construed as medical advice. For any use, the product information guides, inserts, and operation manuals of the drugs and devices should be consulted. Edwards Lifesciences, and their respective affiliates disclaim any liability arising directly or indirectly from the use of drugs, devices, techniques, or procedures described in this reference card.

WARNING: Any reference to X-ray exposure, intravenous contrast dosage, and other medication is intended as a reference guideline only. The guidelines in this document do not substitute for the judgment of a health care provider. Each scan requires medical judgment by the health care provider about exposing the patient to ionizing radiation. Use the As Low As Reasonably Achievable (ALARA) radiation dose principle to balance factors such as the patient's condition, size, and age; region to be imaged; and diagnostic task.

NOTE: Algorithms/protocols included in this reference card are for educational reference only. The authors do not endorse or support any one specific algorithm/protocol. It is up to each individual clinician and institution to select the treatment that is most appropriate.

Philipp Blanke, MD is a paid consultant for Edwards Lifesciences

Edwards, Edwards Lifesciences, and the stylized E logo are trademarks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2021 Edwards Lifesciences Corporation. All rights reserved. DOC-0181405 Rev A

Edwards Lifesciences • One Edwards Way, Irvine, CA 92614 USA • edwards.com

