



# CT Acquisition & Reconstruction for the SAPIEN 3 Valve with Alterra Adaptive Presept




SAPIEN 3 valve



Alterra adaptive presept

# CT Image Acquisition & Reconstruction

Parameters	
Scan coverage	<ul style="list-style-type: none"> <li>• Minimum coverage: RVOT, MPA and branch PAs</li> <li>• Coverage can be extended to include RV and LV if CT is used for chamber quantification (higher radiation exposure!)</li> </ul>
Contrast application	<ul style="list-style-type: none"> <li>• Optimized for right heart/MPA assessment</li> <li>• Biphasic injection protocol using diluted contrast (70%/30% contrast/saline mixture) to reduce streak artifacts, 4-5ml/sec flow rate</li> <li>• Bolus tracking for acquisition timing, ROI placed in RV or pulmonary artery</li> </ul>
Data acquisition mode	<ul style="list-style-type: none"> <li>• Volume CT scanner (GE Revolution, Canon/Toshiba AQ1): ECG-gated volume acquisition</li> <li>• All other scanner: retrospective ECG-gating</li> </ul>
Heart cycle coverage	<ul style="list-style-type: none"> <li>• Systole and diastole should be covered (ideally the entire cardiac cycle), to allow for dynamic assessment of landing zone/potentially contractile RVOT; if dose modulation is required, a peak dose window covering 20%-80% of the RR interval is recommended</li> </ul>
Image reconstruction	<ul style="list-style-type: none"> <li>• Multiphasic reconstruction covering entire cardiac cycle in 5-10% or 50msec increments</li> <li>• 512-512 matrix, field of view limited to cardiac structures, MPA and branch PAs</li> <li>• Maximize use of iterative reconstruction</li> <li>• <math>\leq 0.75\text{mm}</math> slice thickness</li> </ul>




## Systole and Diastole should be covered!

This allows for assessment of:

- Dynamic changes throughout the cardiac cycle, with increase in MPA diameter and length in systole
- RVOT contractility


Largest MPA dimensions throughout cardiac cycle determine anatomical suitability



Scan range needs to cover at least MPA/branch PAs and RVOT.

Avoid truncating these structures!

RV/LV can be included into scan range if need for chamber quantification.



This handout is provided as an educational resource to medical personnel by Edwards Lifesciences (the “Author”). The information in this handout has been compiled from then-currently available literature. Although every effort has been made to faithfully report the information and keep it up to date, the Author cannot be held responsible for the completeness or accuracy. This handout is not intended to be and should not be construed as medical advice. For any use, the product information guides, inserts, and operation manuals of the drugs and devices should be consulted. Edwards Lifesciences, and their respective affiliates disclaim any liability arising directly or indirectly from the use of drugs, devices, techniques, or procedures described in this handout.

**WARNING:** *Any reference to X-ray exposure, intravenous contrast dosage, and other medication is intended as a reference guideline only. The guidelines in this document do not substitute for the judgment of a health care provider. Each scan requires medical judgment by the health care provider about exposing the patient to ionizing radiation. Use the As Low As Reasonably Achievable (ALARA) radiation dose principle to balance factors such as the patient’s condition, size, and age; region to be imaged; and diagnostic task.*

**NOTE:** *Algorithms/protocols included in this handout are for educational reference only. The authors do not endorse or support any one specific algorithm/protocol. It is up to each individual clinician and institution to select the treatment that is most appropriate.*

**CAUTION: Federal (United States) law restricts these devices to sale by or on the order of a physician.**

Edwards, Edwards Lifesciences, the stylized E logo, Alterra, Edwards SAPIEN, Edwards SAPIEN 3, Edwards SAPIEN 3 Ultra, eSheath, SAPIEN, SAPIEN 3, and SAPIEN 3 Ultra are trademarks of Edwards Lifesciences Corporation or its affiliates. All other trademarks or service marks are the property of their respective owners.

© 2022 Edwards Lifesciences Corporation. All rights reserved. DOC-0188837 Rev B

Edwards Lifesciences • edwards.com

One Edwards Way, Irvine CA 92614 USA • +1.949.250.2500 • +1.800.424.3278



Edwards