

Optimising echocardiography

to support aortic stenosis
assessment and severity grading



Edwards

The ongoing challenge of distinguishing severe from less than severe aortic stenosis

- Patients are graded as having severe aortic stenosis (sAS) when they meet both the haemodynamic criteria ($V_{\max} \geq 4$ m/s or mean gradient ≥ 40 mm Hg) and aortic valve area (AVA) criteria (≤ 1 cm²)¹
- However, evidence suggests patients with discordant indices of V_{\max} 3.0–3.9 m/s and AVA ≤ 1 cm² may have sAS if certain criteria apply^{1,2}
 - For these patients, clinical echocardiography (echo) reports characterising the discordant indices as less than severe, could have adverse clinical consequences²

In a retrospective analysis of echocardiographic data from 807 patients with AS²



Patients with discordant indices were found to be less commonly characterised as having sAS on an echocardiography report despite meeting echocardiographic criteria for sAS²



It is vital all patients with sAS are identified to ensure they receive the treatment they need^{3,4}

View subgroups of indices that could indicate sAS



Echocardiography is key for assessing and grading AS

An integrated, stepwise approach that uses multiple indices* is recommended when assessing and grading AS in clinical practice¹



Echocardiography can be used to²

- Confirm the presence of AS
- Assess the degree of valve calcification, left-ventricular (LV) function and wall thickness
- Detect the presence of other associated valve disease or aortic pathology and provides prognostic information



Optimal imaging is needed to determine timing of intervention for sAS, which – once symptomatic – is crucial for survival

*Transvalvular velocity/gradient, AVA, valve morphology, flow rate, LV morphology and function, blood pressure and symptoms.

References: 1. Baumgartner H, et al. Eur Heart J Cardiovasc Imaging. 2017;18(3):254–275. 2. Baumgartner H, et al. Eur Heart J. 2017 Sep 21;38(36):2739–2791.

Optimising assessments with echocardiography



Edwards

A stepwise approach is recommended when grading AS

- Both the EACVI and ESC/EACTS recommend a stepwise approach to assessing AS severity using echocardiography^{1,2}
- As high gradient generally indicates sAS,¹ this guide focusses on the low gradient track and the measurements needed to confirm true sAS:

Step 1:

Calculating mean gradient and peak velocity

Step 2:

Measuring left-ventricular outflow tract (LVOT) and LVOT velocity time integral (VTI) to calculate AVA

Step 3:


Excluding measurement errors to ensure the identified discordant indices are correct i.e., low velocity/gradient where AVA <1cm²

Step 4:

Defining flow status by calculating the stroke volume index (SVi)

Step 5:


Assessing LV ejection fraction (LVEF) to identify patients with true sAS (including the use of dobutamine echocardiography [Step 6] and calcium score [Step 7])



Before taking measurements: Confirming visual signs of AS

Check if the valve has visual signs of AS (i.e. signs of calcification and little exertion of the leaflets) to identify who has sAS and requires further assessment¹

Valve Anatomy



Identify number of cusps in systole, raphe if present

Assess the mobility of cusp and commissural fusion

Assess calcification of valve


Adapted from Ring L, et al. 2021.²

View the visual signs of AS on aortic valve anatomy at different AS severities



Step 1: Measuring mean gradient and peak velocity

To measure mean gradient, obtain the highest velocity using continuous wave (CW) doppler and use multiple acoustic windows (e.g. apical, suprasternal, right-parasternal)



Adapted from Baumgartner H, et al. 2017¹

To improve data acquisition

- Optimise the signal by decreasing gain, increasing wall filter and adjusting baseline curve and scale
- Use grey scale spectral display with expanded time scale
- Adjust the velocity range and baseline so that the velocity signal fits but fills the vertical scale

To optimise measurement

- Measure the maximum velocity at the peak of dense velocity curve. Be careful to avoid noise and fine linear signals
- Trace around the outer edge of the dense signal to obtain VTI
- Calculate mean gradient from the traced velocity curve
- Report window where maximum velocity is obtained

View jet velocity at different AS severities




Pitfalls to avoid when measuring velocity and mean gradient

Avoiding the following pitfalls during CW doppler acquisition and utilising multiple windows will help to optimise assessments

Pitfall to avoid & incorrect example

❌ AS severity can be underestimated if imaging of the vena contracta (the narrowest and highest velocity portion of the jet) is off-axis¹

- The presence of a faint transaortic flow (red trace) may indicate that the insonation beam is angled away from the vena contracta¹



❌ The direction of the transaortic jet may often be anterior and to the right¹

- The severity of AS in up to 50% of patients could be underestimated if the peak velocity is measured from a single apical window¹



Measurement limitations²


VELOCITY/GRADIENT

- Parallel alignment of the ultrasound beam is required for correct measurement
- Flow dependent
- The accuracy of the pressure gradients is dependent on accurate velocity


Correct measurement & example

✅ Adjust the imaging window to find the most dense and uniform continuous-wave spectral profile¹

Use multiple acoustic windows (e.g. apical, suprasternal and right parasternal) during the CW doppler assessment to minimise errors²



✅ For example, in the example here the right parasternal window showed peak velocities that were 0.5 m/sec higher. Demonstrating why the use of multiple acoustic windows is recommended¹




All images Hahn RT, et al. 2018.¹

Step 2a: Measuring LVOT and LVOT VTI

To calculate AVA and other indices

Take LVOT diameter measurements in the parasternal long axis (PLAX) view¹

LVOT diameter CSA



Take LVOT diameter measurements at the same anatomic level as the velocity recording

Diameter is used to calculate a circular cross-sectional area (CSA)


Measure from the inner edge of the septal endocardium and the anterior mitral leaflet in the mid-systole (blue arrow)

Take the measurement parallel and adjacent to the aortic valve or at the site of the velocity measurement

Adapted from Baumgartner H. et al. 2017.¹

Measure LVOT VTI using pulsed wave (PW) doppler and apical long-axis or 5-chamber view

LVOT VTI – PW – Doppler sample



Report window where maximum velocity and mean gradient are obtained

Trace around the outer edge of dense signal to obtain VTI


Measure the maximum velocity at the peak of dense velocity curve. Be careful to avoid noise and fine linear signals

Pitfalls to avoid when measuring LVOT



Step 2b: Calculating AVA

To assess if a patient with low gradient AS is at risk of sAS, calculate the AVA using the continuity equation.



Baumgartner H, et al. 2017¹

- The continuity equation is based on conservation of mass and assumes that the stroke volume ejected through the LVOT all passes through the AVA¹
- The equation requires the following measurements:¹
 - A_1 = LVOT diameter for calculation of circular cross sectional area
 - V_1 = LVOT (peak V & VTI) recorded with PW doppler
 - V_2 = AV (peak V & VTI) recorded with CW doppler

Step 3: Excluding measurement errors

To identify who has moderate AS and might have sAS with discordant indices

A major limitation of the standard continuity equation is the assumption of a circular outflow tract shape, despite LVOT becoming progressively more elliptical (rather than circular) in many patients¹

Consequently, when LVOT diameter is squared for the calculation of cross-sectional area, it has the potential to be the greatest source of measurement error in the continuity equation¹



LVOT VTI

Can give you cardiac output

Stroke volume = LVOT area x LVOT VTI.
Cardiac output is = Heart rate x stroke volume

Step 4: Defining flow status

An AVA $\leq 1.0 \text{ cm}^2$ suggests but does not confirm sAS, therefore it is important to exclude measurement errors and define flow status using SVi¹

- SVi is calculated as: $\text{SVi} = \text{Stroke volume in mL} / \text{Body surface area in m}^2$
- If SVi is normal ($>35 \text{ mL/m}^2$) sAS is unlikely even if AVA is calculated as $\leq 1.0 \text{ cm}^2$
- As many as 35% of sAS patients may be in a low-flow state ($\text{SVi} \leq 35 \text{ mL/m}^2$) and require careful haemodynamic evaluation to determine if they have sAS²

Step 5: Assessing LV ejection fraction To identify patients with true sAS

For patients with low SVi, LVEF should be assessed to identify those with preserved ejection fraction versus those with reduced ejection fraction

- Patients with lower than expected gradients despite preserved LVEF can lead to an underestimation of severity, which may delay aortic valve replacement³
- Use LVEF to determine the cause of the low-flow state:²
 - In patients with a reduced LVEF ($<40\%$) dobutamine stress echocardiography (**Step 6**) can be used to differentiate between true sAS and pseudo-sAS (where uncertainty remains, assessing the calcium score by MSCT (**Step 7**) can help identify true sAS)¹
 - In patients with preserved LVEF an integrated approach, including the use of MSCT (**Step 7**), should be used to determine true sAS



A high calcium score can confirm sAS while a low score makes sAS highly unlikely

Help avoid misdiagnosis of AS severity in your patients by optimising echocardiogram assessment

- sAS can pose unique challenges that require special considerations during work-up and imaging
- Evidence has shown some patients with sAS and discordant indices are being characterised as having less severe AS in clinical echocardiography (echo) reports, which could have adverse clinical consequences¹
- Both the EACVI and ESC/EACTS recommend an integrated, stepwise approach that uses echocardiography to assess multiple indices* when grading AS severity in clinical practice^{2,3}

*Transvalvular velocity/gradient, AVA, valve morphology, flow rate, left ventricular (LV) morphology and function, blood pressure and symptoms.

References: 1. Raddatz MA, et al. Open Heart. 2020 Aug;7(2):e001331. 2. Baumgartner H, et al. Eur Heart J Cardiovasc Imaging. 2017;18(3):254–275. 3. Baumgartner H, et al. Eur Heart J. 2017 Sep 21;38(36):2739–2791. 4. Holmes DR Jr, et al. Ann Thorac Surg. 2012;93(4):1340–1395.

Thank you

Edwards, Edwards Lifesciences, the stylized E logo are trademarks or service marks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

These materials are sponsored solely by Edwards and were not authorized by any other manufacturer.

© 2021 Edwards Lifesciences Corporation. All rights reserved. PP--EU-2360 v1.0

Edwards Lifesciences • Route de l'Etraz 70, 1260 Nyon, Switzerland • edwards.com




Edwards

Subgroups of indices that could indicate sAS^{1,2}

Subgroup	AVA	Mean pressure gradient	SVi	Ejection fraction	Haemodynamic consequences
High gradient AS	<1 cm ²	>40 mmHg	-	-	sAS assumed irrespective of whether LVEF and flow are normal or reduced
Low-flow, low-gradient AS with reduced ejection fraction	<1 cm ²	<40 mmHg	≤35 mL/m ²	<50%	<p>Low-dose dobutamine echocardiography is recommended in this setting to distinguish true from pseudo sAS:</p> <ul style="list-style-type: none"> This is defined by an increase to an AVA of >1.0 cm² with flow normalisation
Low-flow, low-gradient AS with preserved ejection fraction	<1 cm ²	<40 mmHg	≤35 mL/m ²	≥50%	<p>Exclude measurement errors and other reasons for such echocardiographic findings.</p> <p>The degree of valve calcification by MSCT can help identify true sAS:</p> <ul style="list-style-type: none"> A high calcium score can confirm sAS while a low score makes severe AS highly unlikely
Normal-flow, low gradient AS with preserved ejection fraction	<1 cm ²	<40 mmHg	>35 mL/m ²	≥50%	Generally indicates moderate AS.

Pitfalls to avoid when measuring LVOT¹



Adapted from Hahn RT & Pibarot P. 2017.¹

Visual signs of AS in AVA¹



Normal



Aortic Stenosis



Mild-to-moderate
aortic stenosis




Severe aortic stenosis

Adapted from Otto CM. 2008.¹




If the data does not support the visual signs of AS, check for measurement errors


Visual signs of AS in AVA¹




Normal



Aortic Stenosis
<2.5 m/sec



Mild-to-moderate
aortic stenosis 2.5-4 m/sec



Severe aortic stenosis
>4 m/sec

Adapted from Otto CM. 2008.¹



If the data does not support the visual signs of AS, check for measurement errors